

Sports Medicine Evidence Matters

Research Bulletin

Stryker CrossBlade Arthroscopic Cutters Exhibit Lowest Rate of Particle Generation

Top-Level Summary:

The particle generation of four CrossBlade cutters was compared to four competitive arthroscopic cutters. This test was intended to measure the amount of particulate generated when operating the arthroscopic cutter under known conditions. **All four CrossBlade cutters were found to generate significantly less particulate than all competitive cutters tested.**¹

Methods:

The particle generation from eight arthroscopic cutters was measured:

- Stryker CrossBlade Aggressive Max
- Stryker CrossBlade Smooth Bite
- Stryker CrossBlade Tomcat HC
- Stryker CrossBlade Dual Edge
- Arthrex Bone Cutter
- Linvatec Ultracut
- Dyonics Bone Cutter Platinum
- Dyonics Full Radius

Each arthroscopic cutter was operated under controlled conditions of speed, time, and loading. During the test, deionized water was flushed through the cutter removing any debris being generated. Next the water was pumped through a filter capturing particulate for evaluation. Using automated particle counting equipment (Clemex Technologies Inc.) the percent area per million (PAPM) of particles was measured. Significance was assumed at $p \leq 0.05$.

Results:




Figure 1

The CrossBlade Aggressive Max, Smooth Bite, Tomcat HC, and Dual Edge arthroscopic cutters all exhibited significantly less particle generation than all other blades evaluated. (Ave. + st.dev.)¹




Figure 2

Images of particulate gathered on 0.8µm filter paper taken using 200X microscope.

Clinical Relevance:

While inherent to the design of metal arthroscopic cutters, metallic microparticles have been shown to elicit an inflammatory cellular response when found at high concentrations.² **The significant reduction in particle generation by the CrossBlade arthroscopic cutters may result in safer arthroscopic surgery.**

References:

1. Stryker DHD12753 and DHD13100
2. Pedowitz et al. "Arthroscopic surgical tools: A source of metal particles and possible joint damage" *Arthroscopy* 29(9): 1559–1565, 2013.

5670 Greenwood Plaza Blvd. Ste. 200
Greenwood Village, CO 80111
t: 201 831 5000
www.stryker.com
www.sportsmedicine.stryker.com

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks; CrossBlade, Stryker. All other trademarks are trademarks of their respective owners or holders.

1000901580 Rev D
Copyright © 2016 Stryker
Printed in USA