

MV3[®] bariatric bed

Safe and dignified bariatric care

MV3 is a safe and durable bariatric bed designed with a patient's dignity and independence in mind. With versatility in mind, MV3's deck litter expands to three widths to accommodate patients of various sizes. What's more, MV3's low bed height, fall prevention technologies, and a Stryker-exclusive multi-pivotal patient helper encourage interdependence and mobility for patient's of size.



Maximum weight capacity: MV3 is durable and can accommodate 1,100 pound safe working load.

3 Expandable widths: Our expandable width has three settings – 36/42/48 inches – and can accommodate patients of various sizes.

Advanced fall prevention technologies: This includes low bed height, integrated bed exit with pause feature, 360 degree bed lights, and Stryker exclusive iBed Awareness.

Dignified design: MV3 helps preserve the dignity of bariatric patients, as it's been designed to look and feel like the rest of the Stryker fleet.

Features

1. iBed Awareness*
2. Single-zone bed exit alarm
3. One-touch siderails with four positions
4. Universal brake at foot and head end
5. Patient control lock out and total bed motion master lock out
6. Litter deck expansion handle
7. Manual and electronic CPR
8. Digital head of bed
9. Auxiliary cord outlet
10. 360-degree lights on footend corners



*iBed Awareness monitors local bed status information and alerts caregivers visually or audibly if preset parameters are compromised.

Specifications

MV3 Specifications

Model number	5900
Weight capacity	
Safe working load	1102 lbs (500 kg)
Maximum patient weight	1027 lbs. (465.8 kg)
Overall length	92" (233.7 cm)
Patient surface length	80" (203.2 cm)
Width(s) (inside of siderails)	
36" (91.4 cm), 42" (106.7 cm), 48" (121.9 cm)	
Width (s) (outside of siderails)	
40.5" (102.9 cm), 46.5" (118.1 cm), 52.5" (133.4 cm)	
Height range (to litter top)	
High	28" (71.1 cm)
Low (standard with 6" casters)	12" (30.5 cm)
Litter positioning	
Backrest	0° to 60°
Knee gatch	0° to 20°
Trend./reverse Trend.	+12° to -12°
Caster diameter	6" (15 cm)

Bariatric Support Surface

expandable foam	288505550001
-----------------	--------------

MV3 standard features

- Enclosed central locking brake with neutral and steer function
- Brake pedals at head and foot end
- 6" dual casters (helps prevent floor denting)
- One-touch cardiac chair position
- Digital HOB angle indicator—displayed on siderails and footboard
- Robust siderails capable of 4 positions: Up, half-height, low, and tucked
- 3 expandable bed widths
- One-sided dual width expansion
- Obstruction sensors with motion inhibition—centrally located
- Digital bed status indicator on home screen—displays bed exit, iBed, low height, HOB 30+, brake set
- Integrated scale*
- Single zone-expanding bed exit system
- Bed exit pause feature with automatic re-arm
- Stryker-exclusive iBed Awareness Monitoring System
- Integrated pump holder
- One-touch foot elevate position
- IV pole sockets

MV3 optional accessories

- | | |
|--------------------|------------|
| • 2 stage IV pole | 5900000144 |
| • Patient Helper | 5900000145 |
| • O2 bottle holder | 5900000142 |

- 2 foley bag hooks
- 10 restraint strap loops
- CPR—electronic and manual lever
- Removable and expandable headboard
- Integrated patient siderail controls
- Caregiver siderail and footboard bed controls
- One-touch trend/reverse trend positioning
- Foot end digital trend/reverse trend display
- Locking, removable footboard
- Dimmable 360 degree nightlights with 5 levels of programmable intensity
- Patient exit light
- Auto-contour
- Single motion bed lockout for caregiver/patient (CPR overrides)
- Master motion bed lockout for caregiver/patient (CPR overrides)
- DC actuators
- Two-way integrated siderail speakers

* ± 3% of the total patient weight for patients weighing 225 (102 kg) to 1102 lb (500 kg)

This product is intended for patient ≥ 250 lbs and ≥ 30 BMI

Stryker reserves the right to change specifications without notice.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: iBed, MV3, Stryker. All other trademarks are trademarks of their respective owners or holder.

Mkt Lit-1808 12 JUL 2019 Rev B
Copyright © 2021 Stryker

3800 E. Centre Avenue
Portage, MI 49002 USA
t: 269 329 2100
toll free: 800 327 0770

stryker.com

