

Symmetry[®] Plus

Patient Room Recliner



Versatile recliners for a healing environment

Our stylish Patient Room Recliner maximizes patient and guest comfort, while providing uncompromised durability. Choose from Rocker or Wall-Saver models, with open or upholstered arms available on each model.



Step-lock brakes are easy to access to help improve caregiver efficiency.



Rear push bar helps make transporting a patient easier.



Rolling footrest provides added comfort while in transit.

Backrest stops in a variety of recline positions

Posture-correct lumbar and lateral cushioning

Pull-out foot tray gives patients added comfort while in transit

Rocker or Wall-Saver models available

Extended, wide arm caps facilitate easy entry and exit

Symmetry Plus Patient Room Recliners



Standard features

- Stationary or mobile
- Mobile with 3 inch front rigid/rear swivel locking casters
- Mid-ottoman cushion provides extra support when feet are elevated
- Chair construction includes metal seat frame and sinuous wire springs
- Highly resilient polyurethane foam cushions
- Removable back cushion enables thorough cleaning and service
- Pull-out foot tray gives patients added comfort while in transit
- Rigid push bar for easy maneuverability
- Stationary arms and fixed aides provide stability for drainage catheters or electromechanical mounts (open-arm only)
- Egress assist polyurethane arm caps
- Meets flammability standards of BIFMA F1-1978 — revised 1980 and CAL 117
- Meets ANSI/BIFMA X5.4 2005 per tests conducted in a certified test facility

Optional features

- Open arm or upholstered arm
- Egress assist wood arm caps
- Hinged push bar for maneuverability
- Four-wheel swivel with locking casters
- Choose from a variety of stains and fabrics
- CAL 133 in select fabrics
- Foley Bag hooks (Wall-Saver only)
- 3 inch front/rear swivel locking casters

Specifications

	Rocker recliner	Wall-saver
Model number	3500-000-750	3500-000-730
Overall width	30" (76.2 cm)	30" (76.2 cm)
Overall height	45.5" (115.6 cm)	47.5" (120.7 cm)
Depth		
Overall	41.5" (105.4 cm)	40.5" (102.9 cm)
Seat	19.5" (49.5 cm)	19.5" (49.5 cm)
Recline/leg rest up	70" (177.8 cm)	70" (177.8 cm)
Weight	139 lb. (63 kg)	131 lb. (59.4 kg)
Seat width	23" (58.4 cm)	23" (58.4 cm)
Back height from seat	32.5" (82.6 cm)	33.5" (85.1 cm)
Seat height from floor	20.5" (52.1 cm)	20.5" (52.1 cm)
Arm height from floor	27" (68.6 cm)	27" (68.6 cm)
Safe working load	350 lb. (158.8 kg)	350 lb. (158.8 kg)

Warranty

Limited lifetime warranty

Warrantied 350 lb. (158.8 kg) weight capacity

Rocker recliner

Rear lock lever allows chair to remain stationary

Wall-saver

Reclines away from the wall to save space

Stryker reserves the right to change specifications without notice.

Furniture is manufactured exclusively for Stryker Medical.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, Symmetry. All other trademarks are trademarks of their respective owners or holder.

3800 E. Centre Avenue
Portage, MI 49002 USA
t: 269 329 2100
toll free: 800 327 0770

stryker.com

Mkt Lit-85-02022005 Rev H.0
Copyright © 2018 Stryker
Printed in the U.S.A.