

MX-PRO® R3

x-frame ambulance cot



3800 E. Centre Avenue
Portage, MI 49002 USA
t: 269 329 2100
toll free: 800 327 0770
stryker.com

Standard Features

- Transport at load height capability
- Positive action height adjustment
- Easy-to-use release handle design
- Color-coded controls
- Dampened action during hot drops
- High visibility powder-coated frame
- Lightweight, durable aluminum construction
- Scientifically optimized lift bar design
- Lower lifting bar
- Seven height positions
- Integrated bumper system
- Lift-capable safety bar
- Perforated litter surface
- One-hand release breakaway head section with safety bar
- Floor-mounted safety hook
- One-hand release, fold down side rails
- One-hand release, infinite positioning, pneumatically assisted backrest
- Oversized wheels with sealed caster and wheel bearings
- Reflective labeling
- Sealed bolster mattress
- Shock positioning
- Two lap belts and one four-point shoulder restraint
- Single wheel lock

Optional Features

- Heavy duty two- or three-stage IV poles (patient right or left)
- Base lift bar
- Hard or soft base storage
- High height kit for 33 inches
- Defibrillator platform
- Pocketed or nonpocketed head end storage
- Height limit kit*
- Permanent or removable O₂ bottle holders (head or foot end)
- Pull handle
- Sealed flat mattress
- Premium mattress
- Dual wheel lock
- X-frame guards
- Head extension with pillow
- Pillow only
- Rigid head end storage tray
- Equipment hook
- Backrest storage pouch

Specifications

Model Number	6082
Height¹	
Position 1	13.5 in. (34.3 cm)
Position 2	21 in. (53.3 cm)
Position 3	25.5 in. (64.8 cm)
Position 4	29 in. (73.7 cm)
Position 5	32 in. (81.3 cm)
Position 6	35 in. (88.9 cm)
Position 7	37.5 in. (95.3 cm)
Length	
Standard	80.5 in. (204.5 cm)
Minimum	62 in. (157.5 cm)
Width	23 in. (58.4 cm)
Weight²	81 lb. (37.7 kg)
Wheels	
Diameter	6 in. (15.2 cm)
Width	2 in. (5.1 cm)
Articulation	
Backrest	2° – 73°
Shock Position	+ 14°
Maximum Weight Capacity	650 lb (294.5 kg)
Minimum Operators Required	
Occupied Cot	2
Unoccupied Cot	1
Recommended Fastener System	
Floor Mount	Model 6370, 6377 or 6378
Wall Mount	Model 6371
Recommended Loading Height*	Up To 32 in. (81.3 cm)

Warranty

- One year parts and labor or two years parts
- One year parts and labor or fifteen months parts and labor
- Two year parts only
- Lifetime on all welds*
- * Five year service life

*Cot may be loaded from any height and a height limiting kit (p/n 6060-202) is available to limit the load height of the cot.

¹ Height measured from bottom of mattress at seat section to ground level.

² Cot is weighed without mattress and restraints.

Stryker reserves the right to change specifications without notice.

In-service video included with every order.

The MX-PRO R3 is designed to conform to the Federal Specification for the Star-of-Life Ambulance KKK-A-1822. The yellow and black color scheme is a proprietary trademark of the Stryker Corporation.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: MX-PRO, Stryker. All other trademarks are trademarks of their respective owners or holder.

