

Specifications

Model number	6392
Length	
Overall length	70 inches (178 cm)
Length with mass casualty inclusion	80 inches (203 cm)
Width	
Foot end	9 inches (23 cm)
Head end	19 inches (48 cm)
Weight	
Total weight	62.5 pounds (28.5 kg)
Floor plate assembly	16.5 pounds (7.5 kg)
Transfer assembly	46 pounds (21 kg)
Height	6 inches (15 cm)
Minimum operator required	
Occupied cot	Two
Unoccupied cot	One
Recommended loading height	22 inches to 36 inches (56 cm to 91 cm)

Meets dynamic crash standards with Power-PRO XT (AS/NZS-4535, BS EN-1789 with X-restraints and SAE J3027) with X-restraints and Performance-PRO XT (SAE J3027 with X-restraints).

Performance-LOAD is rated IP26 with inductive charging option

Stryker reserves the right to change specifications without notice.

Specifications are rounded to the nearest whole number. Conversions are calculated before rounding.

The Performance-LOAD cot fastener system is designed to conform to the Federal Specification for the Star-of-Life Ambulance KKK-A-1822.

Patents pending.

Technical Support

Our Technical Support comprises a team of professionals available to help with your Performance-LOAD needs. Contact via phone at 1 800 STRYKER or email at medicaltechnicalsupport@stryker.com

Stryker's ProCare Services

Every day, you count on your medical equipment to perform at its best. With ProCare Services, our people help to ensure your equipment is ready to perform when it's needed and make it easier to get the most from your investment. When an issue arises, we promise that we'll work to solve it — performing repairs quickly and correctly in accordance with the terms of the service agreement.

ProCare isn't just a service program. It's a partnership you can count on to give you one less thing to worry about, and one more reason to feel confident you're doing all you can for your clinicians, staff and patients.

All ProCare agreements provide:

- Stryker-authorized service representative
- Stryker-direct factory parts
- Two-hour callback response time
- Fixed service costs up front
- Increased uptime
- 24-72 hour equipment turnaround time*

*Based on the provisions of the Service Agreement and the location of the product.

Flex Financial Program

Our financial programs provide a range of smart alternatives designed to fit your organization's needs. We offer flexibility beyond a cash purchase with payment structures that can be customized to meet budgetary needs and help to build long-term financial stability. Contact your account manager for more information.

Contact your local sales representative or visit stryker.com

Emergency Care

This document is intended solely for the use of healthcare professionals. A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate Stryker's product offerings. A healthcare professional must always refer to operating instructions for complete directions for use indications, contraindications, warnings, cautions, and potential adverse events, before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area. Specifications subject to change without notice.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Flex Financial, Performance-LOAD, Performance-PRO, Power-LOAD, Power-PRO, Power-PRO XT, ProCare, SMRT, Stryker. All other trademarks are trademarks of their respective owners or holders. The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

Mkt Lit-1241 15 DEC 2015 Rev C
Copyright © 2022 Stryker



Manufactured by:

3800 E. Centre Avenue
Portage, MI 49002 USA
P + 1 269 329 2100
Toll free + 1 800 327 0770
stryker.com

Distributed in Canada by:

Stryker Canada
2 Medicorum Place
Waterdown, Ontario
L8B 1W2
Canada
Toll free 800 668 8323

Performance-LOAD[®]

manual fastener



Performance-LOAD

manual fastener

Our safest, most advanced manual cot fastener.*

During the most critical moments of an emergency, you're focused on the patient—and we're focused on you.

The Performance-LOAD manual fastener helps keep both the patient and caregiver safe by guiding and securing the cot during a call.

It's time to make your job easier, safer and more efficient.*

*When used with Power-PRO 2, Power-PRO XT and Performance-PRO XT cots.

- Meets SAE J3027 dynamic crash test standards for enhanced occupant safety
- Guided loading and unloading support improves operator efficiency by reducing cot drift
- Strategically designed with the ability for future modification to support Stryker's Power-LOAD® fastening system
- Optional inductive charging optimizes SMRT™ battery management
- Mass casualty compliant with the inclusion of wall or floor mounts and appropriate hardware
- Power washable and easily removable for maximum decontamination and lowered potential for biohazards



Head end dock
Secures the head end of the cot into place when in transport position.

Foot end latch
Locks the foot end of the cot into place when in transport position.

Adaptable floor plate and safety hook
Compatible with Stryker's Power-LOAD technologies. The safety hook connects to the cot safety bar during loading and unloading.

Linear transfer system
Guides and centers the cot during loading and unloading.

Removable toolkit
Easy to clean and maintain with a simple removal process.

Inductive charging (option)
Charges the SMRT battery when in transport position (no cable or connectors required).

Positive action foot end release
Disengages the cot.



Safety

- Adheres to dynamic crash testing SAE J3027
- Keeps occupants safe in the event of an ambulance accident



Efficiency

- Guides loading and unloading
- Helps reduce cot drift



Reliability

- Life cycle tests simulate use and over the expected life
- Tested and manufactured to the highest standard

Performance-LOAD cot compatibility

The Performance-LOAD compatibility option is available for the Power-PRO 2, Power-PRO XT, Power-PRO IT and Performance-PRO XT. This system meets dynamic crash test standards, enhancing occupant safety for Power-PRO 2, Power-PRO XT and Performance-PRO XT cots. Optional inductive charging automatically charges Stryker's SMRT batteries.



Power-PRO™ 2 cot



Power-PRO XT™ cot



Performance-PRO™ XT cot



Power-PRO™ IT cot