

## Stair-PRO<sup>®</sup>

### Stair Chair



3800 E. Centre Avenue  
Portage, MI 49002 USA  
t: 269 329 2100  
toll free: 800 327 0770  
ems.stryker.com

#### Standard Features

- High visibility powder-coated frame
- Color-coded controls
- Extendable foot end lift handles
- Locking rear lift handles
- Molded hand grips
- Lightweight, rugged aluminum construction
- Oversized rear wheels with sealed bearings
- Dual wheel locks
- Compact storage size
- Positive action locking mechanism
- Power washable
- Grease-free maintenance

#### Model 6252 Standard Features

- Innovative Stair-TREAD™ system
- Extendable upper control handle
- 4 in. (10 cm) front caster wheels

#### Optional Features

- Two-piece molded ABS seat
- One-piece traditional vinyl seat
- Color-coded patient restraints<sup>4</sup> (plastic buckles)
- Non-absorbent patient restraints<sup>4</sup> (metal buckles)
- Polyester patient restraints<sup>4</sup> (metal buckles)
- Patient head support (model 6252 only)
- O2 bottle holder
- Vinyl head support
- Polypropylene head support
- IV clip
- Foot support

#### Specifications<sup>1</sup>

• Models	6252
• Height	37.5 in. (95.25 cm)
• Width	20.5 in. (52 cm)
• Depth <sup>2</sup>	28 in. (71 cm)
• Folded Depth	8 in. (20 cm)
• Weight	32.8 lb (14.9 kg)
• Maximum Load <sup>3</sup>	500 lb (228 kg)

#### Warranty

- One year parts and labor or two year parts only
- fifteen months parts and labor
- Two year parts and labor
- Two year parts, labor and one year soft goods
- Lifetime on all welds\*

#### Extended warranties available.

\*7-year service life

<sup>1</sup>Dimensions are measured from the outermost edges of the main frame. Specifications are rounded to the nearest whole number. Conversions are calculated before rounding.

<sup>2</sup>Depth dimensions are measured with extendable handles retracted.

<sup>3</sup>Maximum load capacity is total weight distributed in accordance to basic human anatomy. EMTs must consider the weight of the patient, equipment and accessories when determining the total weight on the product.

<sup>4</sup>Patient restraint set includes two cross-chest restraints and one ankle restraint. Cross-chest restraints can also be used in a chest and lap restraint configuration.

Stryker reserves the right to change specifications without notice.

Specifications are rounded to the nearest whole number. Conversions are calculated before rounding.

The yellow and black color scheme is a registered trademark of Stryker Corporation

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Stair-PRO, Stryker**, . All other trademarks are trademarks of their respective owners or holder.