

# Safe care for you and your patients

Help address the risk factors of pressure  
injury and healthcare worker injury



# Are your patients at risk for pressure injuries?

Pressure injuries are a significant health issue and one of the biggest challenges your facility faces on a day-to-day basis. Aside from the high cost of treatment, pressure injuries also have a major impact on your patients' lives and on your hospital's ability to provide appropriate care to patients.<sup>1</sup>

## Too common and costly for patients, families, and the healthcare system

Pressure injuries affect more than  
**2.5M** patients per year<sup>2</sup>




The overall prevalence  
of pressure  
injuries is **9.3%**<sup>3</sup>



About

**60,000**

patients die as a direct  
result of a pressure  
injury each year<sup>4</sup>



Cost to treat pressure injuries  
can range from

**\$20,900 –  
\$151,700**

depending on the stage of injury.<sup>2</sup>



**Pressure injuries  
are commonly  
seen in high-risk  
populations:** elderly  
patients who are very ill;  
critical care patients<sup>1</sup>

# The primary goal—**prevention**

As a nurse, you're often on the **frontline of defense against pressure injuries**. Constant vigilance and communication across units is essential to reduce the rate of pressure injuries in high-risk patient populations.

## **A challenging condition**


**Addressing pressure injuries has always been a challenge, not just for caregivers, but also for the healthcare industry, because the epidemiology of pressure injuries varies by clinical setting.<sup>1</sup>**

The **sacrum and the heels** are the most common site to develop pressure injuries.<sup>5</sup>



## **The fearsome four:**

Address pressure injury risk factors<sup>6</sup>



## **Additional risk factors<sup>7</sup>**

- **Reduced mobility or immobility**
- **Acute illness**
- **Extremes of age**
- **Vascular disease**
- **Level of consciousness**





# Are you **at risk for injury?**

Healthcare professionals incur musculoskeletal (MSK) injuries/nonfatal injuries on the job more often than any other private industry occupation.<sup>8</sup>

Healthcare workers are more likely to get injured on the job than any other occupation—more than construction laborers, firefighters, and police officers.<sup>8</sup>

**The most common tasks that lead to injury are patient:**


- **Lifting**
- **Transferring**
- **Repositioning**<sup>9</sup>

Paying the price	Risks on the rise
<p><b>\$15,800</b></p> <p>Average compensation claim due to patient handling<sup>10</sup></p>	<p> <b>Patient obesity levels are projected to increase</b><sup>13</sup></p>
<p><b>\$37,000</b></p> <p>Average direct cost associated with an occupational back injury of a healthcare provider<sup>11</sup></p>	<p> <b>Average age of nurses has risen</b><sup>14</sup></p>
<p><b>\$27,000 –</b></p> <p><b>\$103,000</b></p> <p>Cost of nurse turnover<sup>12</sup></p>	<p> <b>Aging healthcare workers statistically become more vulnerable to musculoskeletal injuries</b><sup>15</sup></p>
	<p> <b>Estimates project 260,000 unfilled nursing jobs by 2025</b><sup>16</sup></p>

## Who takes care of patients **when healthcare workers are injured on the job?**

Many healthcare workers leave the profession early due to debilitating arm, back, and shoulder injuries. More nurses are worried about getting a back injury than contracting an infectious disease.<sup>17</sup> And for good reason:


- Based on an ANA survey **56% of nurses** have experienced MSK pain that was caused by or made worse at work<sup>17</sup>
- **80%** continued to work despite having MSK pain according to an ANA survey<sup>17</sup>
- **48%** of hospital worker injuries that result in days away from work are caused by overexertion or bodily reaction that includes motions such as lifting, bending, reaching, or slipping without falling<sup>18</sup>



# How is your **hospital addressing safe patient handling?**

There is no such thing as safe manual lifting of patients, regardless of body mechanics.<sup>20</sup>

30+ years of research and experience shows that relying on proper body mechanics or manual lifting techniques alone is not effective to reduce back and other MSK injuries.<sup>21</sup>



Nurses can lift a cumulative weight of up to **1.8 tons** during an 8hr shift.<sup>22</sup>

## Hospitals are investing in equipment, but healthcare workers are still getting injured

**82%** of healthcare workers

who sustained an MSK injury did not use facility lifting equipment.<sup>23</sup>

### Why aren't **nursing staff using lifting equipment?**

The equipment may not be easily accessible. It may be located away from patient care areas, moved to another floor, unavailable while being reprocessed, or, worse yet, locked away in a closet.



Sage

# Multi-Position MATS®

## Mobile Air Transfer System

**Safely reposition and transfer patients before and after tilt\* procedures**

Tilt procedures present unique challenges to the perioperative team, including **staff injury, workflow issues, and patient sliding.** You shouldn't have to compromise healthcare worker safety to ensure proper patient positioning in the OR. That's why we designed Multi-Position MATS to fit seamlessly into your workflow so you can achieve the best possible outcomes.

### Imaging capability

Acceptable for use in MRI,\*\* X-Ray, CT Scan, and ultrasounds<sup>26</sup>

**45%** of nurse injury is caused by overexertion<sup>24</sup>

Patients in **tilt positions** are vulnerable to **complications and injury** due to sliding<sup>25</sup>



### Break line

Helps ensure proper placement of the patient

### Wipeable

MAT surface is easy-to-clean material



## Multiple benefits



### Repositioning on the OR table

Easily adjust patient position for surgery



### Post-procedure transfer

Safely transfer patients using air-assist technology




### Chest strap

Additional securement device to the OR table



### Cradle inflation

Cradles patient with 2-phase inflation and raises patient in one smooth motion



### Arm wrap

Helps reduce risk of nerve damage by protecting the arms<sup>27,28</sup> and keeps arms adducted during surgery<sup>29</sup>

\*Multi-Position MATS is intended to be used on the following tilt positions: Trendelenburg, Reverse, Trendelenburg, Lithotomy, and the Lateral tilt positions.


\*\*Prevalon MATS is MRI safe by rationale. The device is made from all non-metal materials; therefore MR safety testing was not performed. Compatibility tests did not show artifacts. Based on rationale, the MATS is electronically non-conductive and non-magnetic.

All testing done on simulated patients.



### Foam pad

Skin friendly, redistributes pressure<sup>28</sup> to help reduce the risk of patient skin injury, and can be easily removed post-surgery<sup>29</sup>



### M² Microclimate Body Pad

Included in the system, it effectively absorbs and locks in moisture to help protect patient's skin while allowing air to flow through



### Rail straps

Keep system secure to the OR table<sup>27</sup>



### HEPA

Equipped filter



### Prevalon Air Pump

Easily fits into the Prevalon Air Pump Cart for transport or in cabinets and shelves in patient room

HEPA filter is tested and certified to perform to HEPA standards and is **99.97% efficient in processing 0.3 µm particles.**<sup>30</sup>

## System secures patient in the following tilt positions:



Lithotomy Trendelenburg



Trendelenburg




Reverse Trendelenburg



Lateral tilt



See it in action



### Multi-Position MATS Mobile Air Transfer System

(1) Mat with foam pad,  
(1) Chest strap

(1) M² Microclimate  
Body Pad

**5 systems/case**

**Reorder #3232**

41 in x 51 in  
Weight limit: 400lb / 181 kg



### M² Microclimate Body Pad

36 in x 51 in

**30 pads/case** (6 bags of 5)

**Reorder #7250**

Sage

# Prevalon MATS®

## Mobile Air Transfer System

### Safely and easily transfer patients

The Prevalon Mobile Air Transfer System (MATS) is designed to provide nurses with a safe and easy way to laterally transfer patients. The system uses a cushion of air to move patients laterally from one surface to another with significantly less pulling and without the need for lifting. It is designed to remain with the patient and provide transfer assistance throughout a hospital stay.



### Imaging capability


Acceptable for use in MRI\*, X-Ray, CT Scan, and ultrasounds<sup>28</sup>



Cradles patient with 2-phase inflation and raises patient in one smooth motion

### Expandable hose

Extends to 10 ft to easily accommodate care settings



Helps protect hose from environmental contamination



Integrated into workflow to improve efficiency by allowing caregiver to focus on patient at the bedside



### Prevalon Air Pump

Easily fits into the Prevalon Air Pump Cart for transport or in cabinets and shelves in patient room

HEPA filter is tested and certified to perform to HEPA standards and is **99.97% efficient in processing 0.3 μm particles.**<sup>30</sup>

\*Prevalon MATS is MR safe by rationale. The device is made from all non-metal materials; therefore MR safety testing was not performed. Compatibility tests did not show artifacts. Based on rationale, the MAT is electronically non-conductive and non-magnetic.



### M<sup>2</sup> Microclimate Body Pad

Effectively absorbs and locks in moisture to protect patient's skin while allowing air to flow through



### Integrated head support

Inflation provides comfortable support for the head and neck

### Easy grab handles

Positioned along outer edges

### Quick connect valve

Provides an easy, secure connection and a quick release


## Sage HalfMATS Mobile Air Transfer System

Easily transfer patients for  
OR specialty procedures




### Lithotomy positioning

- Gynecology
- Urology
- Cystoscopy



### Orthopedic positioning

- Hip fracture
- Total hip replacement



### Prevalon MATS

- (1) Mobile Air Transfer Mat
- (1) M<sup>2</sup> Microclimate Body Pad

**10 systems/case**  
**Reorder #3242**

### Recycled

Ask your sales representative for information about our **Save Simply** program

**Reorder #3242-R**

**Reorder #3244**

TAA Compliant

**Reorder #3247**

Without M<sup>2</sup> Microclimate Body Pad



### HalfMATS

- (1) Mobile Air Transfer Mat
- (1) Microclimate Management Pad

**10 systems/case**  
**Reorder #3230**

39 in x 49 in

Weight limit: 500 lb / 227 kg



### M<sup>2</sup> Microclimate Body Pad

36 in x 51 in

**30 pads/case (6 bags of 5)**  
**Reorder #7250**



### Microclimate Management Pad

23 in x 36 in

**30 pads/case (6 bags of 5)**  
**Reorder #7550**

39 in x 81 in

Weight limit: 1000lb / 454 kg

Sage

# AirTAP LC<sup>®</sup>

## Lift Compatible Patient Repositioning System

### Promoting safety and mobility

We understand the importance of early mobility. To achieve your mobility protocols and comply with guidelines, you need a readily available solution that stays under your patient and fits seamlessly into your workflow. That's why we combined our innovative air assist technology with lift compatibility to promote safety and mobility.



### 3-in-1 mobility

#### 1 In-bed positioning

Helps address risk factors of pressure injury

#### 2 Lateral transfer

Utilizes air-assisted technology to safely transfer patients

#### 3 Vertical transfer

Enables caregiver to safely transfer patient to the bedside chair

#### M<sup>2</sup> Microclimate Body Pad

Effectively absorbs and locks in moisture to protect patient's skin while allowing air to flow through

#### Hi-tech fabric

Promotes a healthy microclimate for the skin and is compatible with low air loss technology

Easy to clean surface



**Hose Protection Sleeve**

Helps protect hose from environmental contamination



**Point of Care Power Switch**

Integrated into workflow to improve efficiency by allowing caregiver to focus on patient at the bedside

#### HEPA

Equipped filter



#### Prevalon Air Pump

Easily fits into the Prevalon Air Pump Cart for transport or in cabinets and shelves in patient room

HEPA filter is tested and certified to perform to HEPA standards and is **99.97% efficient in processing 0.3 μm particles.**<sup>30</sup>



**Integrated head support**

Helps ensure patient's head maintains neutral position while inflated



**30° Body and Anchor Wedge System**

Helps sustain a 30° side-lying position and reduces pressure by offloading sacrum

Minimize manual activity as wedges slide into a turn, reducing the need for a log roll

Minimizes patient migration down the bed

Reduces the need for boosting, minimizing shear and friction



**Easy roll handles**

Encourage proper posture for healthcare worker

**Black positioning handles**

Encourage proper ergonomic position for the healthcare worker when preparing to reposition patient



**Intuitive strap system**



Enables multiple lifting positions

**Quick Connect Valve**

Provides an easy, secure connection and a quick release



**Sage AirTAP LC Lift Compatible Patient Repositioning System**

- (1) Lift Compatible Glide Sheet
- (1) 30° Body and Anchor Wedge System
- (1) M<sup>2</sup> Microclimate Body Pad

**5 systems/case**  
**Reorder #7287**

37 in x 65 in  
Weight limit: 600 lb / 272 kg

**Sage AirTAP LC Lift Compatible Patient Repositioning System Glide Sheet**

- (1) Lift Compatible Glide Sheet
- (1) M<sup>2</sup> Microclimate Body Pad

**10 systems/case**  
**Reorder #7289**



**30° Body and Anchor Wedge System**  
5 pair/case  
**Reorder #7295**



**M<sup>2</sup> Microclimate Body Pad**  
36 in x 51 in  
**30 pads/case** (6 bags of 5)  
**Reorder #7250**

Sage

# Prevalon AirTAP<sup>®</sup>

## Patient Repositioning System

### Safe positioning for healthcare workers and patients

Boosting and repositioning patients in bed is a high-frequency task. The Prevalon AirTAP System uses air-assisted technology to significantly reduce the force required to move patients. To reduce the risk of pressure ulcers, the Prevalon AirTAP System maintains an appropriate microclimate and delivers proper sacral offloading.



#### Hi-tech fabric

Promotes a healthy microclimate for the skin and is compatible with low air loss technology

Easy to clean surface

#### M<sup>2</sup> Microclimate Body Pad

Effectively absorbs and locks in moisture to protect patient's skin while allowing air to flow through



#### Hose Protection Sleeve

Helps protect hose from environmental contamination



#### Point of care power switch

Integrated into workflow to improve efficiency by allowing caregiver to focus on patient at the bedside



#### HEPA Equipped filter

#### Prevalon Air Pump

Easily fits into the Prevalon Air Pump Cart for transport or in cabinets and shelves in patient room

HEPA filter is tested and certified to perform to HEPA standards and is **99.97% efficient in processing 0.3 μm particles.**<sup>30</sup>





**Anchor Wedge System**

- Helps maintain natural side lying position and reduces pressure by offloading sacrum
- Minimizes patient migration down the bed
- Reduces the need for boosting, minimizing shear and friction



**Tapered design**

Helps ensure patient's head maintains neutral position while inflated



**Black positioning handles**

Encourage proper ergonomic position for the healthcare worker when preparing to reposition patient




**Easy roll handles**

Encourage proper posture for healthcare worker

**Quick connect valve**

Provides an easy, secure connection and a quick release



**Prevalon AirTAP Patient Repositioning System**

- |  |   |
|--|---|
| (1) Prevalon AirTAP Glide Sheet          | (1) Prevalon AirTAP Glide Sheet           |
| (1) M <sup>2</sup> Microclimate Body Pad | (6) M <sup>2</sup> Microclimate Body Pads |
| (2) 30° Body Wedges                      | (2) 30° Body Wedges                       |

**5 systems/case**  
**Reorder #7217**  
**TAA #7227**

**3 systems/case**  
**Reorder #7218**  
**TAA #7228**

**Prevalon AirTAP Glide Sheet**

**10 each/case**  
**Reorder #7219**  
**TAA #7229**

36 in x 55 in  
 Weight limit: 550lb / 250 kg



**30° Body and Anchor Wedge System**  
**5 pair/case**  
**Reorder #7295**



**M<sup>2</sup> Microclimate Body Pad**  
 36 in x 51 in  
**30 pads/case** (6 bags of 5)  
**Reorder #7250**

Sage

# Prevalon AirTAP XXL

## Patient Repositioning System

**Reduce risk factors for skin injury and patient handling injury when turning and repositioning larger patients**

Boosting and repositioning patients in bed is a high-frequency task. The Prevalon AirTAP XXL System helps reduce the risk of injury to nurses and other caregivers by providing an ergonomically friendly method of turning and repositioning larger patients.



### Hose Protection Sleeve

Helps protect hose from environmental contamination



### Point of care power switch

Integrated into workflow to improve efficiency by allowing caregiver to focus on patient at the bedside

### Expandable hose

Extends to 10 ft to easily accommodate care settings

### HEPA Equipped filter



### Prevalon Air Pump

Easily fits into the Prevalon Air Pump Cart for transport or in cabinets and shelves in patient room

HEPA filter is tested and certified to perform to HEPA standards and is 99.97% efficient in processing 0.3  $\mu\text{m}$  particles.<sup>30</sup>



### Head support

Inflation provides comfortable support for the head and neck



### 30° Body Wedges

Help maintain natural side lying position and reduce pressure by offloading sacrum

Eliminate the need to log roll a patient to place wedges



### Low-Friction Glide Sheet

Full-length glide sheet provides integrated head and foot support

M<sup>2</sup> Microclimate Body Pad accommodates larger patients

45 inches wide to help accommodate larger patients



### Body wedge tail

Easy glide fabric enables at least two healthcare workers to turn larger patients



### Easy roll handles

Encourage proper posture for healthcare worker

Extend width of glide sheet by 8 inches

### Quick connect valve

Provides an easy, secure connection and a quick release



### Prevalon AirTAP XXL Patient Repositioning System

- (1) Prevalon AirTAP XXL Patient Repositioning Glide Sheet
- (2) Prevalon XXL 30° Body Wedges
- (1) M<sup>2</sup> Microclimate Body Pad

1 system/case  
Reorder #7231

### Prevalon AirTAP XXL Patient Repositioning System Glide Sheet

- (1) Prevalon AirTAP XXL Patient Repositioning Glide Sheet
- (1) M<sup>2</sup> Microclimate Body Pad

5 systems/case  
Reorder #7294



### Prevalon AirTAP XXL 30° Body Wedge

1 pair/case  
Reorder #7293



### M<sup>2</sup> Microclimate Body Pad XXL

45 in x 57 in  
20 pads/case (4 bags of 5)  
Reorder #7260

43" x 81"  
Weight limit: 1200lb / 544 kg

Sage

# Prevalon TAP XL/XXL

## Turn and Position System

### Turning and repositioning of larger patients

This system helps reduce the risk of injury to nurses and other caregivers by providing an ergonomically friendly method of turning and repositioning larger patients.

#### M<sup>2</sup> Microclimate Body Pad


Effectively absorbs and locks in moisture to protect patient's skin while allowing air to flow through



30° Body Wedges


Larger size wedges redistribute pressure for bariatric patients

Velcro strips attach to low-friction glide sheet, keeping body wedges in place under patient



Mattress Cover

Secures to most bariatric hospital beds  
Can be used in place of a fitted/flat hospital sheet to help reduce friction



Low-Friction Glide Sheet

Larger size glide sheet and M<sup>2</sup> Microclimate Body Pad accommodate bariatric patients



#### Prevalon TAP XL

- (1) Mattress Cover
- (1) Low-Friction Glide Sheet with Anti-Shear Strap
- (6) M<sup>2</sup> Microclimate Body Pads
- (2) 30° Body Wedges

**1 system/case**  
**Reorder #7220**

41 in x 51 in  
Weight limit: 800 lb / 362 kg

#### Prevalon TAP XXL

- (1) Mattress Cover
- (1) Low-Friction Glide Sheet with Anti-Shear Strap
- (6) M<sup>2</sup> Microclimate Body Pads
- (2) 30° Body Wedges

**1 system/case**  
**Reorder #7230**

46" x 51"  
Weight limit: 800 lb / 362 kg



#### M<sup>2</sup> Microclimate Body Pad

**XL** 41 in x 57 in  
**20 pads/case**  
(4 bags of 5)  
**Reorder #7255**

**XXL** 45 in x 57 in  
**20 pads/case**  
(4 bags of 5)  
**Reorder #7260**



Sage

# Prevalon TAP 2.0

## Turn and Position System

### Turn and position patients safely

The Prevalon Turn and Position System 2.0, unlike lift slings and plastic slide sheets, stays under the patient at all times. It's always ready to assist with turning, repositioning, and boosting the patient. This makes it possible for nurses and staff to achieve compliance to q2 turning protocol while allowing them to focus on providing the best care and minimizing stress on the patient.

### M<sup>2</sup> Microclimate Body Pad

Effectively absorbs and locks in moisture to protect patient's skin while allowing air to flow through



**Anchor Wedge System**

Helps maintain natural side lying position and reduces pressure by offloading sacrum

Minimizes patient migration down the bed

Reduces the need for boosting, minimizing shear and friction



**Low-Friction Glide Sheet**

Works with anchor wedge system to reduce friction and shear

Dermasuede material grips the microclimate body pad to keep it in place



**Black positioning handles**

Assist with quick, gentle microturns to position the patient at the appropriate angle



**Orange boost straps**

Promote proper body mechanics and reduce reliance on grip strength



36 in x 55 in  
Weight limit: 550 lb / 250 kg

### Prevalon TAP 2.0

- (2) 30° Body Wedges with Anchor Wedge System
- (1) Low-Friction Glide Sheet
- (1) M<sup>2</sup> Microclimate Body Pad

5 systems/case

**Reorder #7201**  
**Reorder #7201-WBS**  
(Without boost straps)

### Low-Friction Glide Sheet

20 each/case **Reorder #7296**

- (2) 30° Body Wedges with Anchor Wedge System
- (1) Low-Friction Glide Sheet
- (6) M<sup>2</sup> Microclimate Body Pads

3 systems/case

**Reorder #7206**  
**Reorder #7206-WBS**  
(Without boost straps)



### 30° Body and Anchor Wedge System

5 pair/case  
**Reorder #7295**



### M<sup>2</sup> Microclimate Body Pad

36 in x 51 in  
30 pads/case (6 bags of 5)  
**Reorder #7250**

Sage

# Prevalon TAP 1.0

## Turn and Position System

### Reduce risk factors for skin injury and patient handling injury

The Prevalon Turn and Position System 1.0 is designed to assist and maintain proper patient positioning to offload the sacrum. It also effectively manages moisture and body heat to protect patients' skin. In addition, it helps reduce the risk of injuries to nurses and caregivers by providing a more ergonomically friendly method of repositioning patients.

#### 30° Body Wedges

Soft foam wedges allow staff to easily position patients for consistent protection and pressure redistribution



Positioning handles

Assist with quick, gentle microturns to position the patient at the appropriate angle

#### Anti-shear strap

Anchors the system to the bed to reduce shearing

#### M<sup>2</sup> Microclimate Body Pad

Effectively absorbs and locks in moisture to protect patient's skin while allowing air to flow through

#### Low Friction Glide Sheet

Works with the 30° Body Wedge System to reduce friction and shear

Dermasuede material grips the microclimate body pad to keep it in place



#### Prevalon TAP 1.0

- (2) 30° Body Wedges
- (1) Low-Friction Glide Sheet
- (1) M<sup>2</sup> Microclimate Body Pad

5 systems/case

Reorder #7200

- (2) 30° Body Wedges
- (1) Low-Friction Glide Sheet
- (6) M<sup>2</sup> Microclimate Body Pads

3 systems/case

Reorder #7205

#### Low-Friction Glide Sheet

20 each/case

Reorder #7298



#### 30° Body Wedges

6 pairs/case

Reorder #7297



#### M<sup>2</sup> Microclimate Body Pad

36" x 51"

30 pads/case (6 bags of 5)

Reorder #7250

36 in x 55 in  
Weight limit: 550 lb / 250 kg

Sage

# Prevalon SPS™

## Seated Positioning System

### Designed to promote early patient mobility and reduce risk factors for patient handling injury

Boosting and repositioning patients in the bedside chair can put clinicians at risk for injury. The Prevalon Seated Positioning System provides an easy option for clinicians to safely glide patients to an optimal upright-seated position without lifting. It is uniquely engineered to keep the seated patient in place, minimizing the need for repetitive boosting and repositioning.



#### One-way glide

Unique material resists forward movement, ensuring patients remain in an optimal position



#### Multi-chamber air cushion

Redistributes pressure and provides comfort and security while allowing patient to shift in the chair

#### Microclimate Management Pad

Contains and absorbs fluid to protect patient's skin while allowing air to flow through



#### Multi-grip handles

Promote proper ergonomics for healthcare worker posture and body mechanics



22 in x 33 in  
Weight limit: 350 lb / 159 kg

#### Prevalon Seated Positioning System

- (1) Seated Positioning Cushion with Fastener Strips
- (1) Microclimate Management Pad

5 systems/case  
Reorder #7530

High Elevation (Above 3000 ft)

5 systems/case  
Reorder #7531



#### Microclimate Management Pad

23 in x 36 in

30 pads/case (6 bags of 5)  
Reorder #7550

Sage

# Prevalon Heel Protectors

The heel is the second most common site for pressure injury development.<sup>6</sup> The Prevalon Heel Protector was specifically designed to help reduce the risk of heel pressure injury while keeping the foot and leg in a neutral position.

## Sage Prevalon Heel Protector



- Offloads the heel
- Sequential Compression Device (SCD) compatible

### Low-friction outer shell

Helps maintain patient's freedom of movement by easily sliding over bed sheets

### SCD compatible

### Dermasuede fabric interior

Gently grips limb so it remains fully offloaded, even when patient is moving

### Closure straps

Help secure foot



## Sage Prevalon Heel Protector



- Offloads the heel
- Reduces plantar flexion contracture risk
- Sequential Compression Device (SCD) compatible

### Contracture strap

Helps prevent plantar flexion contracture

### Rip-stop nylon

Helps maintain patient's freedom of movement by easily sliding over bed sheets

### SCD compatible

### Dermasuede fabric interior

Gently grips limb so it remains fully offloaded, even when patient is moving

### Expandable straps

Stretch to accommodate lower limb edema

No sharp edges or irritating surfaces





- Offloads the heel
- Reduces plantar flexion contracture risk
- Helps prevent lateral rotation, reducing risk of peroneal nerve damage
- Sequential Compression Device (SCD) compatible



**Contracture strap**

Helps prevent plantar flexion contracture

**Rip-stop nylon**

Helps maintain patient's freedom of movement by easily sliding over bed sheets

**SCD compatible**

**Integrated anti-rotation wedge**

Helps prevent lateral foot and leg rotation, reducing the risk of peroneal nerve damage

**Dermasuede fabric interior**

Gently grips limb so it remains fully offloaded, even when patient is moving

**Expandable straps**

Stretch to accommodate lower limb edema

No sharp edges or irritating surfaces



**Prevalon Heel Protector I**

Recommended for patients with calf circumference of: 10 in – 18 in (25 cm – 46 cm)

**8 packages/case**  
**Reorder #7305**



**Prevalon Heel Protector II**

Recommended for patients with calf circumference of: 10 in – 18 in (25 cm – 46 cm)

**8 packages/case**  
**Reorder #7300**  
**2 packages/case**  
**Reorder #7302**



**Foot and Leg Stabilizer Wedge**

For use with reorder #7300 and #7302

**10 packages/case**  
**Reorder #7350**



**Prevalon Heel Protector III with integrated wedge**

Recommended for patients with calf circumference of: 10 in – 18 in (25 cm – 46 cm)

**8 packages/case**  
**Reorder #7355**



**Prevalon Heel Protector XL with integrated wedge**

Recommended for patients with calf circumference of: 18 in – 24 in (46 cm – 61 cm)

**2 packages/case**  
**Reorder #7382**



**Prevalon Heel Protector Petite with integrated wedge**


Recommended for patients with calf circumference of: 6 in – 10 in (15 cm – 25 cm)

**8 packages/case**  
**Reorder #7310**  
**2 packages/case**  
**Reorder #7312**

Sage  
Prevalon

# Air Pump with HEPA Equipped filter

- ✓ Our entire portfolio of air-assisted devices use the same HEPA-equipped air pump.
- ✓ Prevalon Air Pump HEPA filter has a 5-year lifespan and is not maintenance intensive.<sup>34</sup>
- ✓ Sage provides the only HEPA-certified air pump for the continuum of care in hospitals today.



## What is your filter replacement and maintenance schedule?



**Standard air pump**  
Replace air filter every 6 months<sup>36</sup>


VS.



**Sage HEPA Equipped Air Pump**  
Replace air filter every 5 years<sup>36</sup>

### Filter time replacement

(for a facility with 50 air pumps)



					
<b>Prevalon Air Pump 120V</b> 1/case Reorder #7455	<b>HEPA Equipped Replacement Filter</b> 4/case Reorder #7465	<b>Prevalon Air Pump Cart</b> 1 cart/case Reorder #7475	<b>Air Pump Bed Adapter</b> 1/case Reorder #7485	<b>Hose Protection Sleeve (HPS)</b> 50/case Reorder #7460	<b>25' inflation hose</b> 1/case Reorder #7443

# Our commitment to you and your patients

We believe building partnerships can deliver better outcomes. Our goal is to provide you with the products, training, education, and support you need to provide the best possible care for your patients.

## Available resources include:



Ongoing product and process in-servicing and education



FocusRN—a clinical education program that includes speaker programs and online learning.



CustomerOne tracking and reporting



Pre- and post-intervention assessments




Dedicated Medical Science Liaison to discuss evidence-based information with you and your team

## Helping you **drive** **change at your facility**

### **Our CustomerOne Value Analysis Program measures, analyzes, and reports**

Changing practice involves lots of effort and, above all else, data. You may know a change is needed, but you lack the evidence to gain acceptance. Evaluation is critical, but you may not have resources to gather, analyze, and report on your own. We can help. Our exclusive team of **CustomerOne** professionals is your expert resource for customized measurement and data analysis.

**Let us help validate your success!**





# Simple interventions. Extraordinary outcomes.

We are your partner for proven prevention. Our industry-leading products solve real healthcare problems and are backed by clinical evidence. Our products allow you to deliver essential patient care with confidence by helping to address risk factors that can lead to infections, skin injury, and caregiver injuries.

**We are driven to solve real problems and make healthcare better for you and your patients.**



Healthcare workers are critical to delivering essential care and improving patient outcomes. Keeping them safe should be a priority for every hospital. **The goal of the Journey to Zero® program is to achieve zero harm for healthcare workers** – which helps improve care and makes hospitals safer for everyone. The program includes a range of repositioning and transfer devices along with educational services, assessments, and utilization tracking to help promote a culture of safety.

**We're committed to actively partnering with your facility to reduce the risk of patient handling injuries and improving outcomes.**

**Let's start the journey to zero – together.**

\*Zero healthcare worker injuries and zero doubt messages are not guarantees and are aspirational in nature.

**References:** **1.** Quick Safety: An advisory on safety and quality issues, Preventing Pressure Injuries, The Joint Commission, Issue 25, July 2016. **2.** Berlowitz, D., VanDeusen Lukas, C., Parker, V., et al. (content last reviewed Oct. 2014). Preventing Pressure Ulcers in Hospitals: A Toolkit for Improving Quality of Care. Agency for Healthcare Research and Quality, Rockville, MD. Retrieved from <http://www.ahrq.gov/professionals/systems/hospital/pressureulcer toolkit/index.html>. **3.** VanGilder, C., Lachenbruch, C., Algrim-Boyle, C., et al. (2017 Jan./Feb.). The International Pressure Ulcer Prevalence™ Survey: 2006-2015: A 10-Year Pressure Injury Prevalence and Demographic Trend Analysis by Care Setting. J Wound Ostomy Continence Nurs, 44(1), 20-28. **4.** Health Research & Educational Trust (2017, April). Hospital Acquired Pressure Ulcers/Injuries (HAPU/I): 2017. Chicago, IL: Health Research & Educational Trust. Accessed at <http://www.hret-hilm.org>. **5.** Salcido R, Lee A, Ahn C, Heel Pressure Ulcers: Purple Heel and Deep Tissue Injury, Clinical Management Extra, Advances in Skin & Wound Care 2011;24(4);374-380. **6.** National Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment for Pressure Ulcers: Quick Reference Guide. Emily Haesler (Ed.). Cambridge Media: Osborne Park, Western Australia; 2014. **7.** Clinical Practice Guidelines: the use of pressure-relieving devices (beds, mattresses and overlays) for the prevention of pressure ulcers in primary and secondary care. Royal College of Nursing, January 2005. **8.** U.S. Bureau of Labor Statistics, [www.bls.gov/news.release/osh2.nr0.htm](http://www.bls.gov/news.release/osh2.nr0.htm) accessed May 19, 2017. **9.** American Nurses Association, Safe Patient Handling and Mobility: Understanding the benefits of a comprehensive SPHM program 2015. **10.** AON Risk Solutions. Health Care Workers Compensation Barometer. December 2014. **11.** Ghosh T, Occupational Health and Hazards among Health Care Workers, International Journal of Occupational Safety and Health 2013;3(1);1-4. **12.** Occupational Safety and Health Administration (OSHA). Safe Patient Handling Programs: Effectiveness and Cost Savings. Retrieved from <https://www.osha.gov/Publications/OSHA3279.pdf> Accessed October 15, 2015. **13.** Wang, Y., Beydoun, M.A., Liang, L., Caballero, B. and Kumanyika, S. K. (2008). Will All Americans Become Overweight or Obese? Estimating the Progression and Cost of the US Obesity Epidemic. Obesity, 16:2323-2330. **14.** American Association of Colleges of Nursing. (2017 May 18). Nursing Shortage Fact Sheet. Retrieved from <http://www.aacnursing.org/News-Information/Fact-Sheets/Nursing-Shortage>. **15.** Occupational Safety and Health Administration (OSHA). Worker Safety in Your Hospital brochure. **16.** Buerhaus, PI, Auerbach, DI, Staiger, The Recent Surge in Nurse Employment: Causes and Implications, Health Affairs July/August 2009 vol.28 no.4 w657-w668. **17.** American Nurses Association Backgrounder, 2011 ANA Health & Safety Survey: Hazards of the RN work environment. **18.** Occupational Safety and Health Administration (OSHA), Caring for Our Caregivers, Facts About Hospital Worker Safety, September 2013. Available at [https://www.osha.gov/dsg/hospitals/documents/1.2\\_Factbook\\_508.pdf](https://www.osha.gov/dsg/hospitals/documents/1.2_Factbook_508.pdf) Accessed September 26, 2017. **19.** American Nurse Today, Special Report: Preventing Patient-Handling Injuries in Nurses, Vol II (5), May 2016. **20.** Safe Patient Handling and Mobility - American Nursing Association, 2013. Safe Patient Handling and Mobility - American Nursing Association publication, Interprofessional National Standards Across the Care Continuum, Silver Spring, MD 2013. **21.** Occupational Safety and Health Administration (OSHA) Safe Patient Handling: Busting the Myths, Caring for our Caregivers Brochure. Available at [https://www.osha.gov/dsg/hospitals/documents/3.1\\_Mythbusters\\_508.pdf](https://www.osha.gov/dsg/hospitals/documents/3.1_Mythbusters_508.pdf). Accessed October 15, 2015. **22.** American Nurses Association, Position Statement: Safe Patient Handling Movement. Available at <http://nursingworld.org/DocumentVault/GOVA/Federal/Federal-Issues/SPHM.html>. Accessed July 19, 2017. **23.** Centers for Disease Control and Prevention. Occupational Traumatic Injuries Among Workers in Health Care Facilities – United States, 2012-2014. MMWR 2015;64(15);405-410. **24.** U.S. Bureau of Labor Statistics, Monthly Labor Review, Occupational injuries and illnesses among registered nurses, November 2018. **25.** Chitlik A, Safe Positioning for Robotic-Assisted Laparoscopic Prostatectomy, AORN Journal 2011;94(1);37-48. **26.** Data on file, Sage Products LLC. MATS CFG Report-13245-1 v1.0. MRComp TR0428-101, Prevalon MATS Rev 004 Imaging Compatibility Testing Summary, Sage Report P700-16-030R. **27.** Data on file, Sage Products LLC. TP21-0180 Multi-Position MATS Comparative Pressure Mapping. **28.** Data on file, Sage Products LLC. Product Evaluation – TrendeMATS Survey, conducted September 2021. **29.** Data on file, Sage Products LLC. TP21-0163 Multi-Position MATS Patient Migration and Device Securement Justification. **30.** Data on file, Sage Products LLC, IESTRP-CC007.3 testing conducted by APC Filtration Inc. Brentford, ON Canada. June 23-24, 2021. **31.** Bouyer-Ferullo S, Preventing Perioperative Peripheral Nerve Injuries, AORN Journal, 2013 Jan;97(1):110-124. **32.** Data on file, Sage Products LLC. TP21-0010: TrendeMATS Arm Wrap Pressure Map, TP21-0013: TrendeMATS IV Arm Wrap Pressure Map. **33.** Data on file, Sage Products LLC. MATS Exertion Test Report P700-16-025R and MATS Exertion Testing Summary. **34.** Data on file, Sage Products, LLC. TP23-0125. December 2023. **35.** Data on file, Sage Products, LLC. TP19-0006. **36.** Medline Comfort Glide Blower, User Manual.

This document is intended solely for the use of healthcare professionals. A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that healthcare professionals be trained in the use of any particular product before using the product. The information presented is intended to demonstrate Stryker's products. A healthcare professional must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area. Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: AirTAP, CustomerOne, FocusRN, Journey to Zero, Prevalon, PrimaFit, Sage, SPS, Stryker. All other trademarks are trademarks of their respective owners or holders. The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.