


# BIO DBM™

## Fiber Boat



The BIO DBM Fiber Boat, featuring Demineralized Bone Fibers (DBF), is an allograft tissue form comprising 100% human tissue. BIO DBM Fiber Boats are pre-moistened with sterile saline, and are designed to be moldable and malleable, while retaining the ability to soak up bone marrow aspirate (BMA), platelet-rich plasma, blood, or saline.

# BIO DBM Fiber Boat



- **Osteoinductive potential:** Demineralized to help expose natural growth factors.<sup>1,2</sup> Each lot passes osteoinductive testing as release criteria.<sup>1</sup>
- **Osteoconductive:** 100% human bone scaffold with no added carrier<sup>1</sup>
- **Applications:**
  - Posterolateral fusion
  - Posterior cervical fusion
  - Any surgical application that requires a non-structural bone void filler

Part #	Description	Size (LxWxH, mm)
7770050	BIO DBM Fiber Boat, Small	50x25x6
7770100	BIO DBM Fiber Boat, Large	100x25x6

1. Data on file at AlloSource

2. Zhang M, Powers RM, and Wolfenbarger L. Effect(s) of the demineralization process on the osteoinductivity of demineralized bone matrix. J Periodontol. 1997; 68:1085-1092

## Spine division

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: BIO DBM, Stryker. All other trademarks are trademarks of their respective owners or holders. DBF is a registered trademark of TheraCell, Inc.

GEN-SS-12\_28434  
Copyright © 2021 Stryker

Stryker Spine  
2 Pearl Court  
Allendale, NJ 07401  
t: 201 749 8000  
stryker.com