



GraftJacket Now™ Ultra-Thick

Regenerative tissue matrix



GraftJacket Now Ultra-Thick is a pre-hydrated, ready-to-use acellular dermal matrix used for reinforcement, support and protection of soft tissue.

Biomechanical properties

- Offers robust suture retention and material strength^{1,2}
- Graft thickness between 3.0mm – 4.0mm

Biology

- Demonstrated growth factor retention post-processing:
 - Tested for presence of growth factors: bFGF, PDGF-bb, TGF-b and BMP 2.³
- Collagen-rich scaffold to help facilitate tissue regeneration.

Ease of use

- Pre-hydrated to allow for immediate use in operative setting
- Two-year shelf life at ambient temperature



Ready-to-use, 3mm – 4mm thick

Product number	Description
8635UESM	GraftJacket Now™ small, 4x4cm, 3mm - 4mm thick
8635UEMD	GraftJacket Now™ medium, 4x8cm, 3mm-4mm thick
8635UELG	GraftJacket Now™ large, 5x7cm, 3mm-4mm thick

References:

1. AlloSource White Paper M8S0062
2. AlloSource White Paper M8S0100
3. AlloSource White Paper M8S0115

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Stryker.

GraftJacket is a trademark of Allosource.

All other trademarks are trademarks of their respective owners or holders.

1000904048 Rev B

Copyright © 2021 Stryker

**Manufactured by:**

AlloSource
6278 S Troy Cir.
Centennial, CO 80111
t: 800 557 3587

Stryker
Sports Medicine
5670 Greenwood Plaza Blvd.
Ste. 200
Greenwood Village, CO 80111

t: 866 596 2022
www.sportsmedicine.stryker.com