


Connect for greater insight

Smart Equipment Management for Neptune 3 & Neptune S

Level up with data-driven insights

Smart Equipment Management (SEM) is an online, cloud-based application designed to enhance the ownership experience for your Neptune rovers. It provides digital visibility to your data, offering reliability and efficiency for your staff. SEM supports troubleshooting insights, balanced product use, simplified service, and asset tracking.



Connect. Empower. Enhance.


Enabled through your docking station

The SEM application is powered by the docking station and is compatible with the Neptune 3 and Neptune S rovers. The latest data is sent from the rover to the docker at the end of every docking cycle over a wireless internet connection to the cloud.¹



Simple, smart set-up

Initial connection is simple and only requires a wireless internet connection to send data. For security reasons, no patient health information is ever included in the data sent.




Easy access on all your devices

Your data driven insights are displayed within the SEM portal (insightcenter.stryker.com).


Balanced product use

Empowers customers to enhance fleet management by strategically rotating rovers, optimizing utilization and ensuring balanced performance.




Asset tracking

Enables customers with the capability to effortlessly view the assigned fleet, encompassing both rovers and dockers, along with their current locations.



Device maintenance

Provides real-time equipment status and detailed error information to improve awareness and time to resolution.



Maximize uptime, stay connected

At Stryker, we know you value your equipment and want to stay connected to it. That's why we now provide Smart Equipment Management (SEM) for both Neptune 3 and Neptune S with our ProCare Services.¹ You have peace of mind and connecting is simple—with a secure internet connection, you're ready to access actionable data.

Stay connected, stay confident

To learn more, please visit smartequipmentmanagement.stryker.com. To contact your Surgical Technologies sales representative, call 800 253 3210.

ProCare Services

Maximize equipment uptime with ProCare Services. We understand that downtime is costly and impacts your daily operations. Our ProCare Services program is designed to help maximize equipment uptime with preventative maintenance options and repair service in the event of a malfunction. Feel confident knowing that all equipment service meets Stryker's quality repair specifications.

Part number	Product description
0703-001-000W	Neptune 3 Rover ProCare + SEM Prevent - Includes 2 preventative maintenance visits annually
0703-001-002W	Neptune 3 Rover ProCare + SEM Protect - Preventative maintenance visits NOT included
0711-001-000W	Neptune S Rover ProCare + SEM Prevent - Includes 2 preventative maintenance visits annually
0711-001-001W	Neptune S Rover ProCare + SEM Protect - Preventative maintenance visits NOT included

1. Capability is dependent on model number and manufacture date.

Surgical Technologies

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Neptune, Stryker, SEM, and Smart Equipment Management. All other trademarks are trademarks of their respective owners or holders.