

AlloFiber™

Demineralized Bone Fibers



- ▶ **100% cortical fibers**
- ▶ **No carrier**
- ▶ **Moldable handling**

AlloFiber

Demineralized Bone Fibers

AlloFiber is comprised of freeze-dried cortical fibers and can be easily hydrated and mixed with autologous bone, bone marrow aspirate, blood, or saline. It is moldable and has a putty-like handling for packing bone voids and gaps.

► Composition

- 100% cortical bone
- No carrier

► Native growth factors

- Processed to preserve endogenous growth factors, such as BMPs, PDGF-BB, and others

► Handling

- Irrigation resistant
- Interconnected fibers maintain osteoconductive benefits of the graft¹
- Moldable, putty-like handling

Item number	Description
86060010	AlloFiber DBF 1cc
86060025	AlloFiber DBF 2.5cc
86060050	AlloFiber DBF 5cc
86060100	AlloFiber DBF 10cc

Liquid hydration guidelines¹

1cc AlloFiber: 1-2cc Hydration
 2.5cc AlloFiber: 2-3cc Hydration
 5cc AlloFiber: 4-5cc Hydration
 10cc AlloFiber: 8-10cc Hydration

AlloFiber is osteoinductive, osteoconductive, and osteogenic¹
 when hydrated with patient's BMA



Note: Product can be mixed with autologous bone, bone marrow aspirate, blood or saline

References

1. Data on file at MTF

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: AlloFiber, Stryker. All other trademarks are trademarks of their respective owners or holders.

Content ID: AP-015760

Copyright © 2022 Stryker

Processed by:

Musculoskeletal Transplant Foundation
 125 May Street, Edison, NJ 08837
 t: 800.433.6576

stryker.com