



Content & Tray Layout - Wrist Fusion



REF#	Description	SET QTY
A	940090 Wrist Fusion Tray	1
B	940093 Labeling Clip Wrist Fusion Tray	1
1	703920 Handle, Large	2
2	703922 Handle Insert, AO, Ratchet, cannulated	2
3	703896 Drill Bit, AO, Ø2.0mm x 135mm, Scaled	2
4	703691 Drill Bit, AO, Ø2.6mm x 135mm, Scaled	2
5	703663 Screwdriver Blade, AO, T8, Self Retaining	2
6	703880 Screwdriver Blade, AO, T10, Self Retaining	2
7	703898 Tap, AO, For 3.5mm x 122mm Screws	1
8	703899 Tap, AO, For 2.7mm x 102mm Screws	1
9	703900 Tap, AO, For 2.4mm x 102mm Screws	1
10	703684 Drill Guide, 2.0mm Drill, Comp/Polyaxial T8	1

REF#	Description	SET QTY
11	703882 Drill Guide, 2.6mm Drill, Comp/Polyaxial T10	1
12	45-80040 Countersink for T8/T10 screws, AO	1
13	703897 Overdrill, AO, Ø2.7mm x 122mm	1
14	703694 Overdrill, AO, Ø3.5mm x 122mm	1
15	703884 Drill Guide, For 2.7mm independent lag screw T8	1
16	703686 Drill Guide, For 3.5mm independent lag screw T10	1
17	703921 Handle, Medium	1
18	703928 Joystick for T10 screw holes	1
19	703888 SpeedGuide for 2.4/2.7mm Screws, T8 (L = 30mm)	1
20	703886 SpeedGuide for 3.5mm Screws, T10 (L = 30mm)	1
21	703891 SpeedGuide Drill, AO, Ø2.0mm (L = 30mm)	1
22	703894 SpeedGuide Drill, AO, Ø2.6mm (L = 30mm)	1
23	705170 Depth Gauge 2.4mm - 3.5mm VariAx Screws	1
24	703927 Joystick for T8 screw holes	1

Content & Tray Layout - Wrist Fusion Drawer



REF#	Description	SET QTY
1	940094 Wrist Fusion Screw Rack T8	1
2	940095 Wrist Fusion Screw Rack T10	1
3	629579 Wrist Fusion Plate Standard Bend, 120mm	2
4	629589 Wrist Fusion Plate Short Bend, 123mm	2
5	629590 Wrist Fusion Plate Straight, 120mm	1
6	705292 Periosteal elevator with round edge	1
7	703818 Kirschner Wire with Stop (ø2mm / 150mm with ø4.5mm Stop) (Pack of 5)	1
8	390164 K-Wire 1.6 x 150mm (pack of 10)	1
9	703938 Bending Iron for VariAx Plates	2

REF#	Description	Size	Length mm
- 657312 - 657332	Locking Screw	3.5mm	12 - 32
- 657412 - 657432	Non-Locking Screw	3.5mm	12 - 32
- 656310 - 656324	Locking Screws	2.7mm	10 - 24
- 656410 - 656424	Non Locking Screws	2.7mm	10 - 24
- 656010 - 656024	Locking Screw	2.4mm	10 - 24
- 656110 - 656124	Non Locking Screws	2.4mm	10 - 24

REF#	Description	SET QTY
A	1500-0005 Drawer For Modules & Screw Racks	1
B	940096 Labeling Clip Wrist Fusion Drawer	1
C	940092 Wrist Fusion Implant & Instrument Insert	1

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and / or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and / or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: VariAx 2, Stryker. All other trademarks are trademarks of their respective owners or holders.

All products listed above are CE marked.

Content ID: VAX-TLC-6, 08-2015

Copyright © 2015 Stryker



Manufactured by:
Stryker GmbH
Bohnackerweg 1
2545 Selzach
Switzerland
www.stryker.com